

Product Datasheet

Human IL3RA protein (Biotin) (orb759178)

Description	Biotinylated Recombinant Human Interleukin-3 Receptor Subunit Alpha is produced by our Mammalian exp...
Endotoxins	Less than 0.1 ng/ug (1 EU/ug)
Form/Appearance	Lyophilized from a 0.2 um filtered solution of PBS, pH 7.4.
Storage	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
Tag	C-Fc-Avi
Note	For research use only
Uniprot ID	P26951
Expiration Date	6 months from date of receipt.