

## Product Datasheet

### Human PD-L2 protein (Biotin) (orb759183)

<b>Description</b>	Biotinylated Recombinant Human Programmed Cell Death 1 Ligand 2 is produced by our Mammalian express...
<b>Endotoxins</b>	Less than 0.1 ng/ug (1 EU/ug)
<b>Form/Appearance</b>	Lyophilized from a 0.2 um filtered solution of 20mM PB, 150mMNaCl, pH 7.4.
<b>Storage</b>	Lyophilized protein should be stored at -20°C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7°C for 2-7 days. Aliquots of reconstituted samples are stable at -20°C for 3 months.
<b>Tag</b>	C-6His-Avi
<b>Note</b>	For research use only
<b>Uniprot ID</b>	<b>Q9BQ51</b>
<b>Expiration Date</b>	6 months from date of receipt.