

Product Datasheet

Canine PD-L1 / B7-H1 Protein (orb759211)

Description	Canine PD-L1, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Phe 19 - His 238 (Ac...
Reactivity	Canine
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PD-L1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 51.4 kDa. The protein migrates as 65-80 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_001278901.1
Purity	90%
Source	Canine PD-L1, Fc Tag (orb759211) is expressed from human 293 cells (HEK293). It contains AA Phe 19 - His 238 (Accession # E2RKZ5-1).
MW	51.4 kDa
Uniprot ID	E2RKZ5
NCBI	NP_001278901.1
Expiration Date	6 months from date of receipt.