

Product Datasheet

Cynomolgus LIV-1 / SLC39A6 Protein (orb759216)

Description	Cynomolgus LIV-1, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 21 - Ile 30...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	LIV-1
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 34.5 kDa. The protein migrates as 50-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_005586923.1
Purity	95%
Source	Cynomolgus LIV-1, His Tag (orb759216) is expressed from human 293 cells (HEK293). It contains AA Leu 21 - Ile 309 (Accession # XP_005586923.1).
MW	34.5 kDa
NCBI	XP_005586923.1
Expiration Date	6 months from date of receipt.