

Product Datasheet

Free Fatty Acid Microplate Assay Kit (orb759221)

Description Detection and Quantification of Free Fatty Acid Content.

Reactivity All

Range 0.01 mmol/L - 1 mmol/L

Tested Applications FA


Storage Shipped and store at 4°C for 6 months.

Note For research use only

Application notes Detection Method: Colorimetric, Sample Types: Serum, Plasma, Urine, Saliva, Milk, Cell culture medias, Other biological fluids, Platform: Microplate reader 550 nm

Sample Types Serum, Plasma, Urine, Saliva, Milk, Cell culture medias, Other biological fluids

Expiration Date Please enquire.



Representative
Standard
Curve