


Product Datasheet

Mouse IGF-2 ELISA Kit (orb76040)

Description	Mouse IGF-2 ELISA Kit
Reactivity	Mouse
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: A25-E91
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 400-1000 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	20030 MW
Uniprot ID	P09535
Sensitivity	5 pg/ml
Expiration Date	Please enquire.



Mouse IGF-2
PicoKine
ELISA Kit
standard ...