


Product Datasheet

Human IL-27 ELISA Kit (orb76044)

Description	Human IL-27 ELISA Kit
Reactivity	Human
Range	156 pg/ml - 10,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: EB13 (R21-K229) + P28 (F29-P243)
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is 0-400 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	27493 MW
Uniprot ID	Q8NEV9
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human IL-27
PicoKine
ELISA Kit
standard ...