


Product Datasheet

Human CCL8/MCP-2 ELISA Kit (orb76046)

Description	Human CCL8/MCP-2 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: E.coli; Immunogen sequence: Q24-P99
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum is 38-72 pg/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	11246 MW
Uniprot ID	P80075
Sensitivity	1 pg/ml
Expiration Date	Please enquire.



Human CCL8/MCP-2 PicoKine ELISA Kit stan...