


Product Datasheet

Human CCL1/I 309 ELISA Kit (orb76051)

Description	Human CCL1/I 309 ELISA Kit
Reactivity	Human
Range	15.6 pg/ml - 1,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: K24-K96
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	10992 MW
Uniprot ID	P22362
Sensitivity	4 pg/ml
Expiration Date	Please enquire.



Human CCL1
PicoKine
ELISA Kit
standard c...