


Product Datasheet

Human CD26/DPP4 ELISA Kit (orb76055)

Description	Human CD26/DPP4 ELISA Kit
Reactivity	Human
Range	312 pg/ml - 20,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: D34-P766
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:100, concentration in serum and plasma is 3127-734 ng/ml
Sample Types	cell culture supernatants, serum, plasma (heparin, EDTA) , saliva and urine.
MW	88279 MW
Uniprot ID	P27487
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human CD26/DPP4 PicoKine ELISA Kit stand...