


Product Datasheet

Human Decorin / DCN ELISA Kit (orb76060)

Description	Human Decorin / DCN ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: G17-K359
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Application notes	Dilution ratio of 1:1, concentration in serum and plasma is around 2 ng/ml
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA).
MW	39747 MW
Uniprot ID	P07585
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Decorin PicoKine ELISA Kit standar...