

Product Datasheet

Coronavirus antibody (orb761188)

Description	Mouse monoclonal antibody to Coronavirus
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	IF
Immunogen	Attenuated virus
Target	Coronavirus
Preservatives	Purified antibody containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by titration
Expiration Date	12 months from date of receipt.