

Product Datasheet

Lamin A/C antibody (HRP) (orb761198)

Description	Mouse monoclonal antibody against Lamin A/C conjugated to HRP
Species/Host	Mouse
Reactivity	Human, Rat and Mouse
Conjugation	HRP
Tested Applications	WB
Immunogen	Recombinant bacterially expressed human lamin A
Target	Lamin A/C
Preservatives	Purified antibody (from supernatant) containing PBS + 0.01% thimerosal
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P48678 , P02545 , P48679
Dilution Range	Optimal antibody dilution should be determined by titration
Expiration Date	12 months from date of receipt.