

Product Datasheet

RNA polymerase II CTD repeat YSPTSPS (phospho-S5) antibody (orb761200)

Description	Mouse monoclonal antibody to RNA polymerase II CTD repeat YSPTSPS (phospho-S5)
Species/Host	Mouse
Reactivity	A. thaliana, Drosophila, Fungi, Hamster, Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, WB
Immunogen	Synthetic peptide: 10 repeats of YSPTSpPS (Human)
Target	RNA polymerase II CTD repeat YSPTSPS (phospho S5)
Preservatives	Purified antibody in PBS + 0.02% sodium azide
Concentration	0.95mg/ml
Note	For research use only
Application notes	SP2/0 Ag.14
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB: Use at an assay dependent dilution. Predicted molecular weight: 217 kDa
Expiration Date	12 months from date of receipt.