

---

## Product Datasheet

### Hsp105/HSPH1 Antibody (orb76122)

**Description**

Hsp105/HSPH1 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human, Mouse, Rat

**Conjugation**

Unconjugated

**Tested Applications**

FC, ICC, IF, IHC, IHC-Fr, WB

**Immunogen**

A synthetic peptide corresponding to a sequence at the C-terminus of human Hsp105 (713-733aa EVMEWMNNVMNAQAKKSLDQD), different from the related mouse sequence by one amino acid, and different from the related rat sequence by two amino acids.

**Form/Appearance**

Lyophilized

**Concentration**

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

**Storage**

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

**Note**

For research use only

**Application notes**

Tested Species: In-house tested species with positive results. Predicted Species: Species predicted to be fit for the product based on sequence similarities. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

**Isotype**

Rabbit IgG

**Clonality**

Polyclonal

**MW**




96865 MW

**Uniprot ID**
**Q92598**
**Sensitivity**

10pg/ml

**Expiration Date**

12 months from date of receipt.



Western blot analysis of  
Lane 1: Rat  
Ova...