

## Product Datasheet

### Tissue DNA Preparation - Column Kit (orb762373)

<b>Description</b>	Human-Animal-Plant DNA Preparation-Column Kit is designed for rapid and high purity isolation of gen...
<b>Storage</b>	store at ambient temperature
<b>Note</b>	For research use only
<b>Application notes</b>	The Animal-Plant DNA Preparation kit allows extraction of up to 10 µg of highly pure genomic DNA in a size range from 200 bp to 50 kb per preparation. The obtained DNA is suitable for a variety of applications, including real-time PCR, southern blot analysis, genotyping and discovery or validation of SNP/SSR markers.
<b>Expiration Date</b>	12 months from date of receipt.