

## Product Datasheet

### Cynomolgus CD79B Protein (orb762410)

<b>Description</b>	Cynomolgus CD79B, His Tag is expressed from human 293 cells (HEK293). It contains AA Ala 30 - Asp 16...
<b>Reactivity</b>	Monkey
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	CD79B
<b>Preservatives</b>	PBS, 0.2 M Arginine, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 17.2 kDa. The protein migrates as 27-32 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	Ala 30 - Asp 161
<b>Purity</b>	90%
<b>Source</b>	Cynomolgus CD79B, His Tag (orb762410) is expressed from human 293 cells (HEK293). It contains AA Ala 30 - Asp 161 (Accession # A0A2K5WTH8-1).
<b>MW</b>	17.2 kDa
<b>Uniprot ID</b>	<a href="https://www.uniprot.org/entry/A0A2K5WTH8">A0A2K5WTH8</a>
<b>Expiration Date</b>	6 months from date of receipt.