

## Product Datasheet

### Mouse CD79B Protein (orb762412)

<b>Description</b>	Mouse CD79B, His Tag is expressed from human 293 cells (HEK293). It contains AA Val 26 - Asp 158 (Ac...
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	CD79B
<b>Preservatives</b>	PBS, 0.2M Arginine, pH7.4
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-10×His
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 17.0 kDa. The protein migrates as 28-38 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	NP_032365.1
<b>Purity</b>	90%
<b>Source</b>	Mouse CD79B, His Tag (orb762412) is expressed from human 293 cells (HEK293). It contains AA Val 26 - Asp 158 (Accession # P15530-1).
<b>MW</b>	17.0 kDa
<b>Uniprot ID</b>	<a href="#">P15530</a>
<b>NCBI</b>	<a href="#">NP_032365.1</a>
<b>Expiration Date</b>	6 months from date of receipt.