

Product Datasheet

Human Plexin B2 Protein (orb762413)

Description	Human Plexin B2, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 20 - Arg 116...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Plexin B2
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 127.0 kDa. The protein migrates as 150-180 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_036533.2
Purity	90%
Source	Human Plexin B2, His Tag (orb762413) is expressed from human 293 cells (HEK293). It contains AA Leu 20 - Arg 1160 (Accession # O15031-1).
MW	127.0 kDa
Uniprot ID	O15031
NCBI	NP_036533.2
Expiration Date	6 months from date of receipt.