
Product Datasheet

Human Netrin-1 Protein (orb762415)

Description

Human Netrin-1, His Tag is expressed from

Val 22 - Ala 604 ...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Netrin-1

Preservatives

PBS, pH7.4 with NaCl

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 67.5 kDa. The protein migrates as 80-100 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_004813.2

Purity

95%

Source

Human Netrin-1, His Tag (orb762415) is expressed from human 293 cells (HEK293). It contains AA Val 22 - Ala 604 (Accession # O95631-1).

MW

67.5 kDa

Uniprot ID
[O95631](#)
NCBI
[NP_004813.2](#)
Expiration Date

6 months from date of receipt.