

Product Datasheet

Canine TROP-2 / TACSTD2 Protein (orb762416)

Description	Canine TROP-2, His Tag is expressed from human 293 cells (HEK293). It contains AA Gln 39 - Thr 281 (...)
Reactivity	Canine
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TROP-2
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 29.4 kDa. The protein migrates as 40-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_852842.5
Purity	90%
Source	Canine TROP-2, His Tag (orb762416) is expressed from human 293 cells (HEK293). It contains AA Gln 39 - Thr 281 (Accession # XP_852842.5).
MW	29.4 kDa
NCBI	XP_852842.5
Expiration Date	6 months from date of receipt.