

Product Datasheet

Mouse PVRIG Protein (orb762417)

Description	Mouse PVRIG, Mouse IgG2a Fc Tag is expressed from human 293 cells (HEK293). It contains AA Ser 35 - ...
Reactivity	Mouse
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	PVRIG
Preservatives	50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-mFc(mlgG2a)
Note	For research use only
Application notes	This protein carries a mouse IgG2a Fc tag at the C-terminus. The protein has a calculated MW of 41.0 kDa. The protein migrates as 50-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Ser 35 - Asp 165
Purity	95%
Source	Mouse PVRIG, Mouse IgG2a Fc Tag (orb762417) is expressed from human 293 cells (HEK293). It contains AA Ser 35 - Asp 165 (Accession # A0A1B0GS01-1).
MW	41.0 kDa
Uniprot ID	A0A1B0GS01
Expiration Date	6 months from date of receipt.