

Product Datasheet

Cynomolgus PVRIG Protein (orb762418)

Description	Cynomolgus PVRIG, Mouse IgG2a Fc Tag is expressed from human 293 cells (HEK293). It contains AA Thr ...
Reactivity	Monkey
Endotoxins	0.1 EU per µg
Conjugation	Unconjugated
Target	PVRIG
Preservatives	50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-mFc(mlgG2a)
Note	For research use only
Application notes	This protein carries a mouse IgG2a Fc tag at the C-terminus. The protein has a calculated MW of 40.7 kDa. The protein migrates as 45-50 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	XP_005549281.1
Purity	95%
Source	Cynomolgus PVRIG, Mouse IgG2a Fc Tag (orb762418) is expressed from human 293 cells (HEK293). It contains AA Thr 41 - Asp 171 (Accession # A0A2K5WVV8-1).
MW	40.7 kDa
Uniprot ID	A0A2K5WVV8
NCBI	XP_005549281.1
Expiration Date	6 months from date of receipt.