

Product Datasheet

Human TGF-Beta 2 / TGFB2 Protein (orb762420)

Description	Human TGF-Beta 2, Tag Free is expressed from human 293 cells (HEK293). It contains AA Ala 303 - Ser ...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TGF-beta 2
Preservatives	TFA and ACN
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 12.7 kDa. The protein migrates as 13-14 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_003229.1
Purity	95%
Source	Human TGF-Beta 2, Tag Free (orb762420) is expressed from human 293 cells (HEK293). It contains AA Ala 303 - Ser 414 (Accession # P61812-1).
MW	12.7 kDa
Uniprot ID	P61812
NCBI	NP_003229.1
Expiration Date	6 months from date of receipt.