

Product Datasheet

Human TGF-beta 3 / TGFB3 Protein (orb762421)

Description	Human TGF-beta 3, premium grade is expressed from human 293 cells (HEK293). It contains AA Ala 301 -...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TGF-beta 3
Preservatives	50 mM HAC, pH3.2
Form/Appearance	Powder
Storage	-20°C
Tag	Native
Note	For research use only
Application notes	This protein carries no "tag". The protein has a calculated MW of 12.7 kDa. The protein migrates as 12 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_003230.1
Purity	95%
Source	Human TGF-beta 3, Tag Free (orb762421) is expressed from human 293 cells (HEK293). It contains AA Ala 301 - Ser 412 (Accession # P10600-1).
MW	12.7 kDa
Uniprot ID	P10600
NCBI	NP_003230.1
Expiration Date	6 months from date of receipt.