

Product Datasheet

Cynomolgus IL-1 beta / IL-1F2 Protein (orb762424)

Description

Cynomolgus IL-1 beta, His Tag is

	AA Ala 117 - Ser 268 (Acc...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-1 beta
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state; -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower.
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 19.2 kDa. The protein migrates as 20-22 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE).
Purity	95%
Source	Cynomolgus IL-1 beta, His Tag (orb762424) is expressed from E. coli cells. It contains AA Ala 117 - Ser 268 (Accession # P79182-1).
MW	19.2 kDa
Uniprot ID	P79182
Expiration Date	6 months from date of receipt.