
Product Datasheet

Human Flt-3 / Flk-2 Protein (orb762430)

Description

Human Flt-3, His Tag is expressed from

Asn 27 - Ser 543 (Ac...

Reactivity

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

Flt-3

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-10×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 60.4 kDa. The protein migrates as 70-90 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence

NP_004110.2

Purity

95%

Source

Human Flt-3, His Tag (orb762430) is expressed from human 293 cells (HEK293). It contains AA Asn 27 - Ser 543 (Accession # P36888-1).

MW

60.4 kDa

Uniprot ID
[P36888](#)
NCBI
[NP_004110.2](#)
Expiration Date

6 months from date of receipt.