

Product Datasheet

Human MARCO Protein (orb762437)

Description	Human MARCO, His Tag is expressed from human 293 cells (HEK293). It contains AA Met 79 - Val 520 (Ac...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	MARCO
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	N-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the N-terminus. The protein has a calculated MW of 45.9 kDa. The protein migrates as 60-65 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Met 79 - Val 520
Purity	90%
Source	Human MARCO, His Tag (orb762437) is expressed from human 293 cells (HEK293). It contains AA Met 79 - Val 520 (Accession # Q9UEW3-1).
MW	45.9 kDa
Uniprot ID	Q9UEW3
Expiration Date	6 months from date of receipt.