

## Product Datasheet

### Human CD37 Protein (orb762445)

<b>Description</b>	Human CD37, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Arg 112 - Asn 241 (Acc...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	CD37
<b>Preservatives</b>	50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-Fc
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 41.3 kDa. The protein migrates as 50-60 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	NP_001765.1
<b>Purity</b>	95%
<b>Source</b>	Human CD37, Fc Tag (orb762445) is expressed from human 293 cells (HEK293). It contains AA Arg 112 - Asn 241 (Accession # P11049-1).
<b>MW</b>	41.3 kDa
<b>Uniprot ID</b>	<a href="#">P11049</a>
<b>NCBI</b>	<a href="#">NP_001765.1</a>
<b>Expiration Date</b>	6 months from date of receipt.