

Product Datasheet

Human TIGIT Protein (orb762451)

Description	Human TIGIT, His Tag is expressed from human 293 cells (HEK293). It contains AA Met 22 - Pro 141 (Ac...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	TIGIT
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 28.8 kDa. The protein migrates as 30-45 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_776160.2
Purity	90%
Source	Human TIGIT, His Tag (orb762451) is expressed from human 293 cells (HEK293). It contains AA Met 22 - Pro 141 (Accession # Q495A1-1).
MW	28.8 kDa
Uniprot ID	Q495A1
NCBI	NP_776160.2
Expiration Date	6 months from date of receipt.