

Product Datasheet

Human CD200 Quick ELISA Kit (orb762966)

Description	Human CD200 ELISA Kit
Reactivity	Human
Range	31.2 pg/ml - 2,000 pg/ml
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: Q31-G232
Storage	Store at 4°C for 6 months. (Ships with gel ice, can store for up to 3 days in room temperature. Refrigerate upon receipt.)
Note	For research use only
Application notes	90-Minute One-Step Protocol
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	P41217
Sensitivity	10 pg/ml
Expiration Date	Please enquire.

