


Product Datasheet

Human Azurocidin ELISA Kit (orb762969)

Description	Human Azurocidin ELISA Kit
Reactivity	Human
Range	62.5 pg/ml - 4,000 pg/ml
Tested Applications	ELISA
Immunogen	Expression system for standard: NS0; Immunogen sequence: I27-P248
Storage	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
Note	For research use only
Sample Types	cell culture supernatants, cell lysates, serum and plasma (heparin, EDTA)
Uniprot ID	P20160
Sensitivity	10 pg/ml
Expiration Date	Please enquire.



Human Azurocidin PicoKine ELISA Kit stan...