


## Product Datasheet

### Human IFNA2 ELISA Kit (orb763006)

<b>Description</b>	Human IFNA2 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	3.12 pg/ml - 200 pg/ml
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Expression system for standard: E.coli; Immunogen sequence: C24-E189
<b>Storage</b>	Store at 4°C for 6 months, at -20°C for 12 months. Avoid multiple freeze-thaw cycles (Ships with gel ice, can store for up to 3 days in room temperature. Freeze upon receiving.)
<b>Note</b>	For research use only
<b>Application notes</b>	Dilution ratio of 1:1, concentration in serum is less than the lowest standard, 3.1 pg/ml
<b>Sample Types</b>	cell culture supernatants and serum
<b>Uniprot ID</b>	<b>P01563</b>
<b>Sensitivity</b>	1 pg/ml
<b>Expiration Date</b>	Please enquire.



Human IFNA2 PicoKine ELISA Kit standard ...