



Product Datasheet

c-Jun/JUN Antibody (orb763125)


Description	c-Jun/JUN Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence of human c-Jun/JUN (PPHLLPQQMPVQHPR).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	43 kDa
Uniprot ID	P05412
Expiration Date	12 months from date of receipt.



Western blot analysis of c-Jun/JUN using...



IHC analysis of c-Jun/JUN using anti-c-Jun/JUN...



IHC analysis of c-Jun/JUN using anti-c-Jun/JUN...