



Product Datasheet

P2Y12/P2ry12 Antibody (orb763137)

Description	P2Y12/P2ry12 Antibody
Species/Host	Rabbit
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	A synthetic peptide corresponding to a sequence of mouse P2Y12/P2ry12 (NRRPKDKDVTKCSFL).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	40 kDa
Uniprot ID	Q9CPV9
Expiration Date	12 months from date of receipt.



Western blot analysis of P2Y12/P2ry12 us...



Flow Cytometry analysis of HEPA1-6 cells...