



Product Datasheet

EIF4A1 Antibody (monoclonal, 3F11) (orb763208)


Description	EIF4A1 Antibody (monoclonal, 3F11)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	FC, ICC, IF, IHC, WB
Immunogen	A synthetic peptide corresponding to a sequence of human EIF4A1 (LQQIELDLKATQALVLA).
Form/Appearance	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Tested Species: In-house tested species with positive results. Other applications have not been tested. Optimal dilutions should be determined by end users. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.
Isotype	Mouse IgG2a
Clonality	Monoclonal
MW	47 kDa
Uniprot ID	P60842
Expiration Date	12 months from date of receipt.



Western blot analysis of EIF4A1 using an...



IHC analysis of EIF4A1 using anti-EIF4A1...



IHC analysis of EIF4A1 using anti-EIF4A1...