



Product Datasheet

CD15 antibody (FITC) (orb763366)

Description	Mouse monoclonal antibody against CD15 conjugated to FITC
Reactivity	Human
Conjugation	FITC
Tested Applications	FC
Immunogen	WERI-RB-1 retinoblastoma cell line
Target	CD15
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: The reagent is designed for analysis of human blood cells using 4 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (0.4 ml) is sufficient for 100 tests.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with fluorescein isothiocyanate (FITC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Dilution Range	FC4 µl reagent / 100 µl of whole blood or 106 cells in a suspension The content of a vial (04 ml) is sufficient for 100 tests
Expiration Date	12 months from date of receipt.



Separation of human CD15 positive neutro...



Flow cytometry surface staining pattern ...