



Product Datasheet

VSV-G-Tag antibody (orb763376)


Description	Mouse monoclonal antibody to VSV-G-Tag
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB 1:5,000 IP 1:200 IF 1:1,000
Expiration Date	12 months from date of receipt.



IP antibody use: 5ug VSV-G Mouse IgG1 pe...



IF analysis of 293T cells transfected with...



1ug VSV-G fusion protein + Primary antib...