

Product Datasheet

C-Myc-Tag antibody (HRP) (orb763392)

Description	Mouse monoclonal antibody against C-Myc-Tag conjugated to HRP
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	WB 1:5,000-10,000 IP 1:200 IF 1:1,000
Expiration Date	12 months from date of receipt.



Western blot analysis of c-Myc Recombina...