


Product Datasheet

Plant actin antibody (orb763398)

Description	Mouse monoclonal antibody to Plant actin
Species/Host	Mouse
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	43kDa
Dilution Range	WB 1:2,000-5,000
Expiration Date	12 months from date of receipt.



Western blot analysis of Arabidopsis wit...