


Product Datasheet

GAPDH antibody (HRP) (orb763404)

Description	Mouse monoclonal antibody against GAPDH conjugated to HRP
Species/Host	Mouse
Reactivity	Canine, Gallus, Hamster, Human, Insect, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep, Yeast
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	37kDa
Uniprot ID	P04406
Entrez	2597
Dilution Range	WB 1:2,000-5,000 IHC 1:100
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa with orb76...