



Product Datasheet

CA IX antibody (orb763440)


Description	Mouse monoclonal antibody to CA IX
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	35-38kDa
Uniprot ID	Q16790
Entrez	768
Dilution Range	WB 1:3,000 IP 1:200 IHC 1:100-200
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) HeLa, 2) 293...



1: Input: HeLa Cell Lysate 2: IP product...



Immunohistochemical analysis of paraffin...