


## Product Datasheet

### Cox2 antibody (orb763490)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Mouse monoclonal antibody to Cox2   |
| <b>Species/Host</b>        | Mouse   |
| <b>Reactivity</b>          | Human, Mouse, Rat   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | IHC   |
| <b>Immunogen</b>           | Synthetic Peptide   |
| <b>Preservatives</b>       | PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.                                 |
| <b>Concentration</b>       | 1mg/ml  |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG1  |
| <b>Clonality</b>           | Monoclonal  |
| <b>MW</b>                  | 70kDa   |
| <b>Uniprot ID</b>          | <a href="#">P35354</a>  |
| <b>Entrez</b>              | <a href="#">5743</a>  |
| <b>Dilution Range</b>      | IHC 1:100-200   |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |



Immunohistochemical analysis of paraffin...