


Product Datasheet

Rabbit IgG (Fc specific, Heavy Chain) antibody (orb763498)

Description	Mouse monoclonal antibody to Rabbit IgG (Fc specific, Heavy Chain)
Species/Host	Mouse
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	55kDa
Dilution Range	WB 1:5,000-20,000 IHC 1:200-500
Expiration Date	12 months from date of receipt.



Western blot analysis of Reduced Rabbit ...