


Product Datasheet

beta-tubulin antibody (HRP) (orb763505)

Description	Mouse monoclonal antibody against beta-tubulin conjugated to HRP
Species/Host	Mouse
Reactivity	Canine, Gallus, Hamster, Human, Insect, Monkey, Mouse, Rabbit, Rat, Sheep, Yeast
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	55kDa
Uniprot ID	P07437
Entrez	203068
Dilution Range	WB 1:20,000-640,000 IHC 1:1000-1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa, with β -t...