


Product Datasheet

Rubisco (Large Chain) antibody (orb763534)

Description	Mouse monoclonal antibody to Rubisco (Large Chain)
Species/Host	Mouse
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	53kDa
Dilution Range	WB 1:2,000-10,000
Expiration Date	12 months from date of receipt.



Western blot analysis of Arabidopsis with...