



Product Datasheet

CaMKII beta/gamma/delta (phospho-Thr287) antibody (orb763540)

Description	Mouse monoclonal antibody to CaMKII beta/gamma/delta (phospho-Thr287)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	50kDa
Uniprot ID	Q13555 , Q13557 , Q13554
Entrez	817 , 816 , 818
Dilution Range	IHC 1:100-200
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...