



Product Datasheet

Beclin 1 \hat{A} antibody (orb763570)


Description	Mouse monoclonal antibody to Beclin 1 \hat{A}
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Concentration	1mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	60kDa
Uniprot ID	Q14457
Entrez	8678
Dilution Range	WB 1:1,000-2,000 IHC 1:100-200
Expiration Date	12 months from date of receipt.



Western blot analysis of 1) 293T Cell Ly...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...