



Product Datasheet

Flag-Tag antibody (orb763619)

Description	Rabbit polyclonal antibody to Flag-Tag
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic Peptide
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Dilution Range	WB 1:5,000
Expiration Date	12 months from date of receipt.



Transfected 293 cell lysate with Flag-ta...



1ug Flag fusion protein + Primary antibo...