


Product Datasheet

Plant-actin antibody (orb763634)

Description	Rabbit polyclonal antibody to Plant-actin
Species/Host	Rabbit
Reactivity	Plant
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant Protein
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43kDa
Dilution Range	WB 1:3,000
Expiration Date	12 months from date of receipt.



Western blot analysis of Nicotiana tabac...